EN: Indication for use

Oil for covering of medium during in vitro fertilisation and micro-manipulation procedures.

Product Description

OVOIL™ is a sterile, light paraffin oil.

For use after equilibration at +37°C and 5 % CO₂, 6 % CO₂ or ambient atmosphere, depending on the intended use.

Storage instructions and stability

Store dark at +2 to +8°C.

OVOIL™ is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

The following is the general procedure for using OVOIL™. For more information on recommended use of OVOIL™, please see the "G Series Manual™" or the "Vitrolife Fertility Systems Manual". The Manuals are available on request from Vitrolife, and may be ordered through the company's internet website, www.vitrolife.com. Directions for use in other languages can also be found at the website.

Remove OVOIL™ from vial using a sterile, non-toxic pipette. Do not pour.

Using a sterile, non-toxic pipette, add culture medium to a vessel labeled with patient ID. For large volume culture, add 0.5 to 1.0 mL to a sterile tube or well and cover entirely with OVOIL™. For small volume culture, place 10 to 100 µL droplets in a sterile culture dish and cover entirely with OVOIL™.

Specifications

STERILE A	SAL 10 ⁻³
Mouse Embryo Assay	
one-cell, expanded blastocyst on day 5 [%]	≥80
one-cell, cell number on day 5	PASS
Micro droplet Mouse Embryo Assay	
one-cell, expanded blastocyst on day 5 [%]	PASS
one-cell, cell number on day 5	PASS
Human Sperm Oil Assay	
Bacterial endotoxins (LAL assay)	
[EU/mL]	< 0.25

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

Symbols

STERILE A	Sterilized using aseptic processing techniques
↗↘ +37°C	Temperature limitation
⊘	Do not re-use, discard after procedure.
📋	Use by - see label.
⚠	Caution: Consult accompanying documents.
REF	Catalog number
LOT	Batch code

Precautions

Do not use OVOIL™ if it appears cloudy.

OVOIL™ does not contain antibiotics.

Vitrolife recommends that media be opened and used only with aseptic technique.

The risks of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Not for injection.

Caution: Federal (US) law restricts this device to sale by or on the order of a physician.

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Vitrolife Sweden AB, Box 9080 (Gustaf Werners gata 2), SE-400 92 Göteborg, Sweden.
Vitrolife Inc., 3601 South Inca Street, Englewood CO 80110, USA.
Vitrolife Pty. Ltd., 55 Flemington Rd, North Melbourne VIC 3051, Australia.
Manufactured by Vitrolife Sweden AB Göteborg, Sweden